

## Real-time 3D optical imaging

High speed, video-rate 3D imaging with a high spatial resolution like confocal microscopes

### Overview

Various 3D imaging methods including X-ray CT, ultrasonography, and confocal laser scanning microscopy are employed in industrial applications, depending on its usage. However, X-ray CT requires the rotation of a light source or the irradiation area to acquire 3D images, which makes rapid imaging difficult. Ultrasonography enables 3D measurement using 2D array elements but suffers from limited spatial resolution. In addition, confocal laser scanning microscopy requires the sequential acquisition of 2D images while moving the observation plane, which makes 3D observation at video rate impossible.

In recent years, 3D cameras and OCT are known as emerging 3D image acquisition techniques. However, the camera-based technique exploiting patterned- or line-illuminations provides limited spatial resolution, and it cannot visualize the internal information of samples. Besides, OCT is, in principle, applicable to reflected light signals only and impossible to measure fluorescence signals.

A series of inventions solves the above issues and provides an imaging method with a simple implementation that enables the acquisition of precise 3D images in **real-time**.

### Product Application

- ❑ Non-destructive inspection device for high-speed detection of extremely small defect
- ❑ Biological sample observation
- ❑ Observation of fast 3D microstructural deformation, etc.

### IP Data

IP No. : JP6537153, US10816474, DE112017002847.7  
 JP2019-117233, PCT/JP2021/031166  
 Admin No. : T15-198, T17-106, T21-072

Please refer to the following pages for each technology detail

Basic technology  
T15-198




Z to X  
transformation  
T17-106



To be nice for  
CW visible laser  
T21-072

### Performance



### Non-destructive 3D imaging comparison table

	This techno.	X-ray CT*	Ultra-sonic*	3D camera*	OCT*
Drawing speed	Fast 60fps	Slow	Slow	Few 10 fps	Few 10 fps
Intrusion depth	Few $\mu$ m	10 cm	Few cm	No. 0 nm	Few mm
Resolution	Few 100nm	Few cm	Few mm	Few Mm	Few $\mu$ m
Color coding	✓		✓	✓	
Exposed to radiation	No.	Yes.	No.	No.	No.

### Contact

\* as reference

# Optical information apparatus and microscope system

Rapid obtaining of optical information in the depth direction to create 3D optical image!

## Overview

Laser scanning microscopy, used for observing a sample by scanning a laser focal spot and detecting reflecting, scattering, or fluorescence signals from a target can be increased acquisition speed of 2D images by high-speed raster scanning of a laser beam. But the rapid acquisition of 3D images is substantially restricted due to the need of changing the moving observation plane.


Using the present invention can be possible the optical 3D imaging acquired at once without moving the optical system or target, because optical information along the optical axis obtained in different lateral position on the detector by separately concentrating the lights .

## Product Application


- ❑ Rapid 3D image acquisition required domain
- ❑ Biological, organism-related and medical field such as bio-function analysis
- ❑ Industrial fields where the fine functional materials are developed such as metals and chemicals

## IP Data


IP No. : JP6537153, US10816474, DE112017002847.7  
 Inventor : KOZAWA Yuichi  
 Admin No. : T15-198



Fluorescent nano particles



Biological cells



**Rapid acquisition of 3D information from 2D-raster scanning of a needle spot**

- Particle tracking
- Video-rate 3D acquisition
- ❑ Possible to obtain rapidly optical information of the target sample from different positions along the optical axis
- ❑ Possible to create 3D optical image in real-time

## Related Works

- [1] Y. Kozawa and S. Sato, Sci Rep 9, 11687 (2019)  
<https://doi.org/10.1038/s41598-019-48265-3>
- [2] Tohoku Univ. New Technology Presentation Meeting 2020, JST [Seminar Video \(Japanese Speaking\)](#)

## Contact

# Holographic optical element and its manufacturing method

## Conversion of depth information linearly into planar information!

### Overview

In 3D imaging, it is important to effectively encode z information, which is the image formation displacement of the target object, into XY information, which is displacement from the optical axis of the image surface. On the other hand, in single-pixel imaging using only a single detector, it is important to encode in efficient manner the XY information of the target into the time information. However, the conventional encoding method has a limit on distance resolution, imaging distance range, z information decoding accuracy and encryption uniqueness.

This invention provides a new and improved hologram optical element, its manufacturing method and an optical device able to mutually convert Z & XY information with simple configuration. This invention is able to realize high accuracy and high spatial resolution for high-speed 3D imaging (patent US10816474) by converting the object depth direction information of the optical detection system into the planar direction information at the detecting surface.


### Product Application

- Optical field such as optical microscope
- Fast sensing for depth direction in 3D space


### IP Data

IP No. : JP7021772  
 Inventor : NAKAMURA T., IGARASHI S., KOZAWA Y.  
 Admin No. : T17-106


### Concept of linear z-x transformation



### Demonstrated linear shift behavior



### Video-rate 3D acquisition



- The image position in the detecting surface shifts linearly in the planar direction according to the object depth position.
- Axial position is retrieved without further calibration processes.

### Related Works

- [1] T. Nakamura et al., Opt. Lett. **43**, 5949 (2018).
- [2] Y. Kozawa et al., Biomed. Opt. Express **13**, 1702 (2022)
- [3] [Press release \(Tohoku Univ.\)](#)

### Contact

## Spot light generator, optical information detection device and microscope

### Fast 3D imaging without the need of multi-photon excitation

#### Overview

High-speed 3D imaging achieved by scanning a non-diffractive, needle-like spot provides high XYZ resolution. The light-needle-spot is produced by a **Bessel beam** with extended **side lobes**, which is known to cause artifacts in the acquired images.

To reduce the influence by the side lobes, a multiphoton excitation process including two-photon excitation is generally used for fluorescence imaging. Accordingly, a microscope system requires a near-infrared, ultrashort pulsed laser source, which significantly rises in the initial cost. In addition, the need for the multiphoton excitation process limits observation targets to fluorescence imaging only.

This invention provides a method to generate a needle spot using a simple configuration that can **suppress side lobes**, without two-photon excitation, using a **visible light source** available on the market. This invention can be applied to high-speed 3D imaging for surface inspection and profiling by detecting reflected or scattered light signals in addition to the fluorescence signal.


#### Product Application

- Non-destructive inspection to detect & remove refined 3D defects, etc.
- Bioimaging
- Surface profiling

#### IP Data

IP No. : PCT/JP2021/031166  
Inventor : KOZAWA Yuichi  
Admin No. : T21-072

Bessel beam under single-photon excitation



Side lobes of Bessel beams

This invention

Suppression of side lobes for single photon excitation

#### Features

- Cost-cut
- Simple implementation
- Suitable for mass production

Since the application has not yet been published, the information can be disclosed after concluding a fee based contract including a confidentiality clause.

#### Related Works

- [1] Japan Society of Applied Physics 82<sup>nd</sup> autumn meeting 2021 presentation document
- [2] Patent US10816474 → [T15-198](#)

#### Contact

Download OnePager



Contact

<https://www.t-technoarch.co.jp/en/contact.html>



Check Out Our Inventions

<https://www.t-technoarch.co.jp/en/anken.php>



Follow us

<https://www.linkedin.com/company/tohoku-techno-arch>



# Leading you to Successful Industrialization



**TOHOKU TECHNO ARCH**

株式会社 東北テクノアーチ