


# Embolizing agent and kit for blood vessel Embolization

No dosage restrictions and no catheter adhesion



## Overview

Vascular embolization treatment, a minimally invasive surgery for various blood vessel - related conditions, has emerged as a crucial method in treating such as hemorrhage, arteriovenous malformation, aneurysms, and hypervascular tumors. Liquid agents present unique advantages, penetrating distal peripheral vasculature inaccessible to solid counterparts while allowing modulation of penetration depth through compositional adjustments. Conventionally used liquid embolic substances include N-butyl cyanoacrylate (NBCA) and ethylene vinyl alcohol copolymer (Onyx™). NBCA embolizes blood vessels by polymerization reaction with blood, and Onyx™ embolizes blood vessels by precipitation and aggregation in blood. However, in NBCA and Onyx™, there is a problem that the embolic substances may stick to the tip of the catheter, injure the inner wall of the blood vessel and the catheter insertion site, and cause bleeding.

According to the present invention, it is possible to provide an embolic agent and a kit for blood vessel embolization comprising the embolic agent and the embolic agent having excellent injectability with a syringe and safety and stability in vivo while having blood vessel embolization properties equivalent to those of conventional embolic agents. The present invention comprises as an active ingredient a copolymer of a monomer containing a cationic functional group and a monomer containing an aromatic group, wherein at least a part of the copolymer has a structure represented by general formula.



Concentration of embolic substance  $\geq 30$  mg/mL & stable blood gel formed 1 min after addition without adhesion to catheter



## Product Application

- Intravascular embolic substances for cerebral intravascular
- Surgical hemostatic agents, adhesives
- Wound dressing agent for skin defect

## Related Works

[1] Z.Jin et al., PNAS 2022 Vol.119 No.42 e2206685119.

## IP Data

IP No. : WO2023/163170  
 Inventor : Hailong FAN, OSANAI Toshiya, Zhiping JIN, KUROKAWA Takayuki, NONOYAMA Takayuki  
 Admin No. : HK24-005

## Contact



**Tohoku Techno Arch Co., Ltd.**  
 Please visit [CONTACT](#) here

Download OnePager



Contact

<https://www.t-technoarch.co.jp/en/contact.html>



Check Out Our Inventions

<https://www.t-technoarch.co.jp/en/anken.php>



Follow us

<https://www.linkedin.com/company/tohoku-techno-arch>



# Leading you to Successful Industrialization



**TOHOKU TECHNO ARCH**

株式会社 東北テクノアーチ